

## Buffer Preparation (Shi Lab)

### 1. 1 M Tris-HCl Buffers

pH	Volume (L)	TrisBase (g)	HCl (ml)
pH 7.0	2	242.2	150-155
pH 7.5	2	242.2	120-125
pH 8.0	2	242.2	80-85

Autoclavable.

### 2. EDTA 0.5 M (pH8.0)

0.5M, 1L: 148 g EDTA  
+ ~30-40 g NaOH to adjust pH  
(or 186 g EDTA-Na.2H<sub>2</sub>O + ~20 g NaOH)

Note: pH adjusted by NaOH is essential for solubility. Autoclavable.

### 3. TAE DNA Electrophoresis Buffer (50 X)

(2 M Tris, 50 mM EDTA)

4 L  
968 g Tris  
228.4 ml glacial acetic acid  
400 ml 0.5 M EDTA 8.0

To make 1x TAE 20 L, add 400 ml 50X buffer into 19.6 L ddH<sub>2</sub>O.

### 4. SDS-PAGE Gel Solutions

	Vol (L)	Tris (g)	HCl (ml)	10% SDS (ml)
4x Lower gel buffer <i>1.5 M Tris-Cl, pH 8.8, 0.4% SDS</i>	2	363.3	50-60	80 ml
4x Upper gel buffer <i>0.5 M Tris-Cl, pH 6.8, 0.4% SDS</i>	2	121.1	70-80	80 ml

#### 4.1 10% SDS

2L:

200g SDS into 2 L, heat to 68°C for solubility. pH ~6.6.

## 5. 5X SDS Loading Sample Buffer

100 ml

	Stock solution	Add volume
250 mM TrisHCl pH6.8	1 M	25 ml
10% SDS		10 g
30% Glycerol		30 ml
<i>5% <math>\beta</math>-mercapitaethanol (or 0.5M DTT)</i>		<i>5 ml</i>
0.02% bromophenol blue	1%	2 ml

## 6. 6X DNA loading sample buffer:

(40% sucrose, 0.01-0.02% BPB)

100 ml

Add 40 g sucrose to 50 ml ddH<sub>2</sub>O, add 2 ml 1% BPB solution, adjust to 100 ml.

## 7. SDS-PAGE Electrophoresis Running Buffer (10x)

*(1x: 25 mM Tris, 192 mM glycine, 0.1% SDS, pH8.3)*

10 L.

303 g Trisbase (FW 121.1)

1440 g glycine (FW 75.07)

100 g SDS

No need to adjust pH

## 8. Transfer Buffer without SDS (10x)

*(1x: 25 mM Tris, 192 mM glycine, pH8.3)*

*10 L*

*303 g Trisbase,*

*1440 g glycine*

*No need to adjust pH*

## 8.1 Transfer Buffer (1x)

500 ml

50 ml of 10x SDS-PAGE running buffer

100 ml of Methanol (final **20% methanol**)

350 ml ddH<sub>2</sub>O

**9. TBS (10x)**

*(1x: 150 mM NaCl, 10 mM Tris pH8.0)*

10 L  
876.6 g NaCl (FW 58.44),  
121.1 g Tris,  
~50-60 ml HCl  
to pH8.0

**9.1 TBS-T (1x)**

20L  
2L 10x TBS  
200 ml 10% Tween20 (final 0.1% v/v)  
ddH<sub>2</sub>O to 20 L

**9.2 Block buffer**

(5% Nonfat milk in TBS-T)  
5g milk in 100 ml TBST

**10. NaCl 4 M**

2 L: 467.5 g NaCl. Autoclavable.

**11. NaOH 10 M**

0.5 L: 200 g

**12. NaAc 3 M**

500 ml: add 204 g NaAc.3H<sub>2</sub>O (FW 136), adjust pH by glacial acetic acid (~60 ml)  
to pH5.2. Autoclavable.

**13. MgCl<sub>2</sub> 1M**

500 ml: Add 101.65 g MgCl<sub>2</sub>.6H<sub>2</sub>O into 500 ml ddH<sub>2</sub>O. Autoclavable.

**14. CaCl<sub>2</sub> 1M**

400 ml: Add 58.8 g CaCl<sub>2</sub>.2H<sub>2</sub>O (FW 147), filter for sterilization.  
Dilute 10x to make 100 mM CaCl<sub>2</sub>.

**15. MgSO<sub>4</sub> 1M**

500 ml: Add 123.3 g MgSO<sub>4</sub>.7H<sub>2</sub>O into 500 ml ddH<sub>2</sub>O. Autoclavable.

**16. ZnCl 0.1M**

250 ml: 3.4 g ZnCl.

## **Stock in -20°C**

### **1. IPTG (1 M)**

1 g IPTG (FW 238.3) resolved in 4.2 ml (~4 ml) ddH<sub>2</sub>O, filter through 0.22 µm filters, aliquot 1 ml in eppendorf. Store at -20°C.

### **2. DTT (1 M)**

5 g DTT (FW 154.25) resolved in 32.5 ml (~30 ml) 10 mM NaAc (pH 5.2), filter through 0.22 µm filters, aliquot 1 ml in eppendorf. Store at -20°C.

### **3. X-gal (20mg/ml)**

Add 5 ml (~4.8 ml) **DMSO** into 100 mg X-gal bottom (FW 408.24). Store at -20°C.

### **4. PMSF (100 mM, =17.4 mg/ml)**

Resolve 1.74g PMSF (MW 174) in **isopropanol**, total 100 ml. Aliquot and store at -20°C or R.T..

### **5. Carbencillin or Ampicillin (50 mg/ml) in water. 1000x**

2.5 g 50 ml.

### **6. Kanamycin (10 mg/ml) in water. 200x**

0.5 g 50 ml.

### **7. Chloramphenicol (34 mg/ml) in ethanol. 200x**

1.7 g/ 50 ml.

### **8. lysozyme 50 mg/ml, 1000x.**

2.5 g/ 50 ml.

### **9. TSA (MW 303):**

Add 1.32 ml Ethanol into each vial (1 mg?) to make the TSA stock 2.5 mM, 5000x. Final concentration of TSA in the cell culture is 0.5 µM (~150 ng/ml).

## **Solutions.**

### **1. Bacteria lysis buffer (GST pull-down binding buffer)**

(50 mM Tris 7.5, 150 mM NaCl, 0.05% NP-40.)

1L

50 ml 1M Tris HCl 7.5;

37.5 ml 4 M NaCl;

5 ml 10% NP-40.

ddH<sub>2</sub>O to 1L.

#### **1.1. GST pull-down binding buffer (1 M)**

(50 mM Tris 7.5, 300 mM NaCl, 0.05% NP-40.)

1L

50 ml 1M Tris HCl 7.5;

75 ml 4 M NaCl;

5 ml 10% NP-40.

ddH<sub>2</sub>O to 1L.

#### **1.2. GST pull-down binding buffer (1 M)**

(50 mM Tris 7.5, 1 M NaCl, 1% NP-40.)

500 ml

25 ml 1M Tris HCl 7.5;

125 ml 4 M NaCl;

50 ml 10% NP-40.

ddH<sub>2</sub>O to 500ml.

## **2. RIPA Buffer**

(50 mM TrisHCl pH7.4, 150 mM NaCl, 2 mM EDTA, 1% NP-40, 0.1% SDS)

1L

50 ml 1 M Tris 7.4,

37.5 ml 4 M NaCl,

4 ml 0.5 M EDTA,

10 ml NP-40.

10 ml 10% SDS.

## **3. Cell Lysis Buffer (Flag-IP buffer)**

(50 mM TrisHCl pH7.4, 250 mM NaCl, 0.5% Triton X100, 10% glycerol, 1 mM DTT, PMSF, PI (Roche))

1L

50 ml 1 M Tris 7.4,

62.5 ml 4 M NaCl,

5 ml Triton X-100,

1 ml 1 M DTT,

100 ml glycerol.

## ChIP Solutions

<b>ChIP Sweeling buffer</b>	<b>Stock</b>	<b>Vol for 500 ml</b>
5 mM PIPES pH 8.0	0.5 M	5 ml
85 mM KCl	3 M	14 ml
1% NP40	10%	50 ml
ddH2O		433 ml
<b>ChIP Nuclei Lysis buffer</b>		
		<b>500 ml</b>
50 mM Tris-Cl pH 8.0	1 M	25 ml
10 mM EDTA	0.5 M	10 ml
1% SDS	10%	50 ml
ddH2O		425 ml
<b>ChIP Dilution buffer</b>		
		<b>500 ml</b>
16.7 mM Tris-Cl pH 8.0	1 M	8.4 ml
167 mM NaCl	4 M	20.8 ml
0.01% SDS	10%	0.5 ml
1.1% Triton X 100	10%	55 ml
1.2 mM EDTA	0.5 M	1.2 ml
ddH2O		414 ml
<b>ChIP Dialysis buffer -Rabbit</b>		
		<b>1000 ml</b>
50 mM Tris-Cl pH 8.0	1 M	50 ml
0.2% Sarkosyl	20%	10 ml
2 mM EDTA	0.5 M	4 ml
ddH2O		926 ml
<b>ChIP Dialysis buffer -Mouse</b>		
		<b>1000 ml</b>
50 mM Tris-Cl pH 8.0	1 M	50 ml
2 mM EDTA	0.5 M	4 ml
ddH2O		946 ml
<b>ChIP Wash buffer-Rabbit</b>		
		<b>1000 ml</b>
100 mM Tris, pH 9.0	1 M	100 ml
500 mM LiCl (MW 42.4)		21.2 g
1% NP40	10%	100 ml
1% Deoxycholic acid (sodium salt. MW 414.5)		10 g
<b>ChIP Wash buffer-mouse</b>		
		<b>1000 ml</b>
100 mM Tris, pH 8.0,	1 M	100 ml
500 mM LiCl (MW 42.4)		21.2 g
1% NP40	10%	100 ml
1% Deoxycholic acid (sodium salt. MW 414.5)		10 g
<b>ChIP Elution buffer</b>		
		<b>Make fresh</b>
50 mM NaHCO3	1 M	
1% SDS	10%	

1.25 M Glycine 200ml  
(MW=75) 18.8 g

0.5M PIPES 200 ml  
30.2g PIPES, 17-19 ml  
10 M NaOH

1 M NaHCO3 (MW  
84) 4.2g/50 ml

5 ml elution buffer:  
0.5 ml 10% SDS  
21 mg NaHCO3