

Gastrointestinal Stem Cell Culture Protocol
Modified KMS 06 May 2014
HPW Lab

Passage for maintenance (24-well plates)

Thaw Matrigel in ice water bath

Thaw Y-27632 and TGFbeta inhibitors (used at 1:1000)

Make sure have enough thawed 50% CM

Split 1:3 every third day

Aspirate CM

Wash wells with PBS-EDTA

Add ~~200~~ 300 μ l trypsin-EDTA/TrypLE

Scratch and suspend Matrigel (with P1000/P200 pipettes)

Incubate plates at 37°C for 5 min

Add 500-800 μ l washing medium and dissociate organoids by vigorous pipetting in plate

Transfer organoid suspension into a 15 ml tube with 5 ml washing medium pipet mix

Centrifuge at **1000** rpm for 5 min (pellet organoids)

Carefully aspirate sup with ~200ul left and note size of pellets for the next step

Add 500ul-1ml washing medium and re-suspend well

Transfer appropriate volume of suspension to a 1.5 ml tube

Centrifuge at 2-3 min (microcentrifuge, 650xg, 4°C)

Or at 8 min (microcentrifuge, 350xg, 4°C)

Aspirate supernatant as completely as possible (finish using pipet)

Suspend organoids in Matrigel (15ul per well)

Do not dilute Matrigel more than 1:3 or 1:2 or it won't solidify

Keep tubes on ice while suspend all samples for plating

Follow Plating Protocol in Crypt Culture Protocol

Amounts needed:

Matrigel, 15 ul per well

6 wells per sample

4 standard samples (FS, SS, FE, SE)

Therefore need a minimum of 360ul Matrigel in total for plating

Use minimum 90ul per sample, so can suspend each sample in 100ul (to avoid air bubbles)

Add final media to the plate once MG is solidified, for an entire 24 well plate you need:

12.5 mL 50% CM Media

12.5 μ L of Y-27632

12.5 μ L of TGFbeta inhibitor