

Lentivirus Packaging in 96-Well Plate

Plasmids: pCMV-DR8.2 (1µg/µl), pCMV-VSVG (1µg/µl), transfer plasmid (40ng/µl)

Transfection reagent: Lipofectamine 2000

Virus harvesting medium:

DMEM + 10% FBS + 1g/ml BSA + 1% P/S

500 ml Hyclone DMEM high glucose

50 ml FBS (Hyclone SH30071.03)

64 ml BSA (10g/100ml in DMEM, Hyclone BSA, SH3057402)

6ml P/S (antibiotics)

1. Continually pass 293T cells in 1:2 for three times before seeding the cells into 96-well plate;
Day 0

2. Dilute the cells to a final concentration of 5×10^5 cells/ml in antibiotics-free culture medium;

3. Seeding 100µl 293T cells into each well of 96-well plate (i.e., 5×10^4 cells/well);

4. Incubate the cells for overnight (about 16~17 hours) at tissue culture incubator;

Day 1

5. Dilute the transfer plasmid to a final concentration of 25 ng/µl;

6. Transfer 4 µl of each transfer plasmid to the corresponding well in a new 96-well plate;

7. Prepare packing plasmid DNA (pCMV-DR8.2 + pCMV-VSVG) mixture as shown below
Add 11µg (pCMV-DR8.2), 1.1µg (pCMV-VSVG), 2.31 ml OPTI MEM to a 5ml-tube @ mix well @ add 21µl to each well of the above new 96-well plate;

8. Dilute lipofectamine 2000 with OPTI MEM

Add 55 µl lipofectamine 2000 to 2.75 ml OPTI medium @ mix well @ incubate at room temperature for 5 min;

9. Add 25 µl to each well of the new 96-well plate containing plasmids;

10. Incubate at room temperature for 20 min;

11. Transfer 50 µl of DNA-lipofectamine mixture to the cells in each well of 96-well plate;

12. Centrifuge the cells at 2000 rpm (about 890xg) for 30 min at room temperature;

13. Incubate the cells at incubator for overnight;

Day 2

14. Carefully remove the culture medium (~120 µl) from each well after 18 hours culture;

15. Add 200 µl BSA-containing DMEM medium to each well;

16. Incubate the cells at tissue culture incubator for 48 hours;

Day 4

17. Spin down the cells at 2000rpm for 10 min;

18. Carefully transfer 190 µl of culture medium from each well into the corresponding new well of a new 96-well plate;

19. Pool 5 µl from each well together for later virus titer;

20. Seal the plate with sterile sealing film and label the plate;

21. Store the plate at -80 C before cell infection;